


ParaMount ECLS Carrier



The ParaMount allows the market leading CentriMag™ system and a Paragon® PMP Oxygenator to be carried safely & quickly during use - making it the lightest, portable ECMO system on the market.

The ParaMount ECLS Carrier represents a compact and mobile ECMO therapy for the future and, presents a simplified answer to patient movement. Integration with the patient's bed allows equipment and the circuit to move together as a single unit, rather than two separate items.



A complete system is achievable with a CentriMag™ blood pump, CentriMag™ V2 console and a Paragon® PMP Oxygenator.

Innovative and space saving ideas were applied to all aspects of the device. It's construction of aluminium and foam encapsulated carbon fibre ensures it delivers an optimum combination of being light weight, compact and strong. This provides an ideal solution to the potentially awkward inter-department patient movement while on circulatory support.