


Adult Tubing Pack with Paragon Adult Maxi PMP Oxygenator



Chalice has over a decade and a half of experience in producing premium tubing packs for extracorporeal cardiopulmonary bypass circulatory support.


Using state of the art manufacturing equipment in our validated class-7 cleanroom production rooms; Chalice only use components and tubing of the highest quality, which is thoroughly checked and supplied by regulated providers from around the world.



Our experience in producing premium tubing packs allows us to acknowledge and adhere to the stringent guidelines set out by the industry's regulatory bodies, while simultaneously satisfying our customers' highest standards and requirements.

The Adult Tubing Pack with Paragon Adult Maxi PMP Oxygenator is supplied as a convenient singly packaged system and contains:

- A coated Paragon® Adult Maxi PMP Oxygenator, with blood inlet and outlet lines pre-connected,
- A 1000ml Soft Shell Reservoir bag with quick release connectors for easy removal after priming,
- Color coded 3/8" arterial and venous tubing lines contained in a sterile barrier AV Loop Tray,
- A pre-connected Oxygen gas line with integrated 0.2 μ gas line filter,



Adult Tubing Pack with Paragon Adult Maxi PMP Oxygenator: Product Code: CME40004